


AngE™DIABETIC

Comprehensive Vascular Screening.

- ✓ Simultaneous Mesurement in under 3 Minutes
- Determine oABI, TBI, Toe Pressure, Pulse Wave Velocity, Heart Rate Variability and more
- ✓ PWI[™] Pulse Wave Index
- ✓ Passive and active Stress Tests





6-Channel Simultaneous Pulse Wave Recording on Wrists, Ankles and Toes using the TOPP-Method (Tissue Optical Perfusion Pressure).



Advanced Vascular Diagnostics

www.sot-medical.com



Comprehensive Screening

The innovative method provides a fast and secure diagnosis of the peripheral vessel status of diabetics.

It allows measuring on ankle and wrists as well as on fingers and toes, including forefoot- and toe pressure measurement.

AngE DIABETIC combines four measuring cuffs and two optical sensors to enable the measurement on **six measuring spots** simultaneously. This creates diverse possibilities to detect blood flow disorders within the terminal vessels as well as to assess the wound healing success of diabetic patients.

Comprehensive Report


The one-page report combines pulse wave and measurement parameters on a single page at the touch of a button.

Simple 3-Minute Measurement

The measurement can be conducted easily by following only a few steps and without significant stress for the patient.

After applying the cuffs on ankles and wrists, as well as the optical sensors on the toes, the system applies a pressure of 180mmHg and decreases it by steps of 10mmHg.

By recording the optical sensors, the examiner can immediately determine the pressure step at which the patient's toes show the first pulsations. Different key indicators, like the oABI, the PWI, the amplitude or the peak time of the pulse wave, are recorded simultaneously.



"The high sensitivity of the optical sensors allow for a good documentation of the pulse waves, even with marginal blood flow. Given the virtually unfiltered display of pulse curves, dicrotic waves can be clearly identified for healthy and elastic arteries."

Dr. Alfred Obermayer

Head of Institute of Functional Phlebologic Surgery, Karl Landsteiner Society



See how it's applied

Simply scan with your smartphone camera and open the link to the video.



